RAW MILK, PASTEURIZED DAIRY PRODUCTS AND CONTAINERS

Sample Collection:

Date: 04/11/2011	Time:	14:00	Collector:	Kennet	h Sykes	6
Temperature Controls:	Raw:	37 °F	Processed:	36 °F	Size:	Quart
Processor/Distributor:	MAOL	A MILK &	ICE CREAM		ID#:	37-43

Sample Receipt:

Date: 04/12/2011	Time: 08:15		
Temperature:	Raw: 1.0 °C	Processed:	1.0 °C
Received by:	Darneice Lyons		

Environmental Microbiology

Sample Group: ES041211-0148

	SAMPLE INFORMATION			RAW AND PROCESSED DAIRY PRODUCTS					CONTAIN		
Lab No.	Code Date	Size	Description	Residual Phos	Inhibitor	Petrifilm Coliform Count	Petrifilm Aerobic Count	DMSCC	PRBC	PRCC	
547	4-28	Half Pint	Skim (paper)	Not Found	Not Found	<1 EPCC/ml	<250 EPAC/ml				
548	4-26	Half Pint	Skim Chocolate (plastic)	Not Found	Not Found	<1 EHSCC/mL	<250 EPAC/ml				
549	4-28	Half Pint	1% (paper)	Not Found	Not Found	<1 EPCC/ml	<250 EPAC/ml				
550	4-28	Half Pint	1% Chocolate (paper)	Not Found	Not Found	<1 EHSCC/mL	<250 EPAC/ml				
551	4-26	Half Pint	1% Strawberry (paper)	Not Found		<1 EHSCC/mL	<250 EPAC/ml				
552	4-28	14 Ounces	2%	Not Found	Not Found	<1 EPCC/ml	<250 EPAC/ml				
553	4-28	14 Ounces	Whole (Homo)	Not Found	Not Found	<1 EPCC/ml	<250 EPAC/ml				
554	4-28	1/2 Gallon	Half & Half	Not Found	Not Found	<1 EHSCC/g	<250 EPAC/g				
555	4-28	1/2 Gallon	Whipping Cream	Not Found	Not Found	<1 EHSCC/g	<250 EPAC/g				
556	5-6	Quart	NF Buttermilk	Not Found		<1 EHSCC/g					
557	5-9	1/2 Gallon	Buttermilk Blend	Not Found		<1 EHSCC/g					
558	4-11	2 Ounces	Past Cream (PT8-39*)	Not Found	Not Found	13 HSCC/g	<250 EPAC/g				
559	4-11	2 Ounces	In Plant Raw (RT1-35*)		Not Found	-	3800 PAC/ml				
560	4-11	2 Ounces	In Plant Raw (RT2-39*)		Not Found		6500 PAC/ml				
561	4-11	2 Ounces	In Plant Raw (RT3-37*)		Not Found		5800 PAC/ml				

Analyzed by:	Darneice Lyons		
Plating Date:	04/12/2011	Time:	09:30
Temperature when analyzed:		1.0 °C	;
SSF:			

Approved By: Susan Beasley

Firean Baaley

CONTROLS:

PAC: 0 Air Density: 0 /15 min	Equip: 0 Diluent and NB:	0/-			
Incubation temperature:	32.0 °C	32.0 °C			
Inhibitor test used:	Delvo P5				
Inhibitor Positive Control:	Purple				
Inhibitor Negative Control:	Yellow				