

Sample Submission and Reporting Form

RAW MILK, PASTEURIZED DAIRY PRODUCTS AND CONTAINERS

Sample Collection:

Date: 09/26/2011 **Time:** 14:15 **Collector:** Ken Sykes
Temperature Controls: **Raw:** 34 °F **Processed:** 36 °F **Size:** Quart
Processor/Distributor: MAOLA MILK & ICE CREAM **ID#:** 37-43

Sample Receipt:

Date: 09/27/2011 **Time:** 08:00
Temperature: **Raw:** 0.5 °C **Processed:** 0.5 °C
Received by: Susan Beasley

Environmental Microbiology

Sample Group: ES092811-0159

SAMPLE INFORMATION				RAW AND PROCESSED DAIRY PRODUCTS				CONTAINERS	
Lab No.	Code Date	Size	Description	Residual Phos	Inhibitor	Petrifilm Coliform Count	Petrifilm Aerobic Count	DMSCC	PRBC PRCC
1377	10-13	Half Pint	Skim (plastic)	Not Found	Not Found	<1 EPCC/ml	360 EPAC/ml		
1378	10-13	Half Pint	Skim Chocolate (plastic)	Not Found	Not Found	<1 EHSCC/mL	580 PAC/ml		
1379	10-10	Half Pint	1% (paper)	Not Found	Not Found	<1 EPCC/ml	260 PAC/ml		
1380	10-8	Half Pint	1% Strawberry (paper)	Not Found		<1 EHSCC/mL	930 PAC/ml		
1381	10-9	Half Pint	1% Vanilla (paper)	Not Found		<1 EHSCC/mL	<250 EPAC/ml		
1382	10-13	Half Pint	2% (paper)	Not Found	Not Found	<1 EPCC/ml	<2500 EPAC/ml		
1383	10-13	14 Ounces	Homo	Not Found	Not Found	<1 EPCC/ml	880 PAC/ml		
1384	10-13	1/2 Gallon	Half & Half	Not Found	Not Found	<1 EHSCC/g	800 PAC/g		
1385	10-13	1/2 Gallon	Whipping Cream	Not Found	Not Found	<1 EHSCC/g	350 PAC/g		
1386	10-23	Quart	Nonfat Buttermilk	Not Found		<1 EHSCC/g			
1387	10-23	1/2 Gallon	Buttermilk Blend	Not Found		<1 EHSCC/g			
1388	9-26	2 Ounces	Past. Cream (PT 8 36)	Not Found	Not Found	<1 EHSCC/g	<250 EPAC/g		
1389	9-26	2 Ounces	In Plant Raw (RT 1 35)		Not Found		21000 PAC/ml		
1390	9-26	2 Ounces	In Plant Raw (RT 2 36)		Not Found		15000 PAC/ml		
1391	9-26	2 Ounces	In Plant Raw (RT 3 34)		Not Found		17000 PAC/ml		

CONTROLS:

Analyzed by: Susan Beasley
Plating Date: 09/27/2011 **Time:** 10:05
Temperature when analyzed: 0.5 °C
SSF:

PAC: 0 **Equip:** 1
Air Density: 3 /15 min **Diluent and NB:** 0 / -
Incubation temperature: 32.0 °C
Inhibitor test used: Delvo P5
Inhibitor Positive Control: Purple
Inhibitor Negative Control: Yellow

Approved By: Susan Beasley